



Product catalogue 2019/2020



# This is TePe GOOD™



For more than 50 years we've worked closely with dental experts to deliver the best dental products available. And we've been successful. Now we're taking the next step towards a healthier life. We're launching TePe GOOD™ – a new sustainable product range and an overall sustainability initiative to make a difference for the good of the planet and for your good.



## GOOD sources

Bio-based raw materials, such as sugar cane, castor oil and wood fibre.



## Green energy

All our products are produced using green energy from renewable resources, partly from our own rooftop solar panel facility.



## Sustainability quality

We are dedicated to improving people's quality of life and take environmental responsibility in all aspects of our work.



## [www.tepe.com/GOOD](http://www.tepe.com/GOOD)

Visit our site for more information about TePe GOOD™, the products, our initiatives and much more.

# Introducing TePe GOOD™ collection.

Introducing TePe GOOD™ – a unique collection of products with the same TePe quality, design and feeling as always, but made from bio-based plastic. Using good, renewable raw materials, we can recirculate 95% of the CO<sub>2</sub> emissions during their life cycle. Naturally, the products are produced using green energy, partly from our own rooftop solar panel facility.

## TePe GOOD™ products



**TePe GOOD™ Toothbrushes**  
Regular Soft, Compact Soft,  
Mini Extra Soft

page 19



**TePe GOOD**  
Mini Flosser™

page 15

**TePe GOOD™**  
Tongue Cleaner

page 27

Bio-based  
plastic

95%  
recycled  
CO<sub>2</sub>

For more information please visit [tepe.com/good](https://tepe.com/good)

## Healthy teeth throughout life

Based on the vision of healthy teeth for all throughout life, TePe works for raised oral hygiene awareness and prevention of oral disease. Our wide product range is developed to make good oral health possible for everyone.



## Interdental Cleaning

---

TePe® Interdental Brush Original.....	8
TePe® Interdental Brush Extra soft.....	9
TePe Angle™.....	10
Accessories.....	12
Practice Box.....	12
TePe® Gingival Gel.....	14
TePe® Dental Tape.....	14
TePe GOOD Mini Flosser™.....	15
TePe EasyPick™.....	16

## Toothbrushing

---

TePe GOOD™ Regular.....	19
TePe GOOD™ Compact.....	19
TePe GOOD™ Mini.....	19
TePe Select™.....	20
TePe Select™ Compact.....	20
TePe Mini™.....	20
TePe Supreme™.....	20
TePe Supreme™ Compact.....	20
TePe Nova™.....	20
TePe Graphic™.....	21
TePe Kids™.....	21
TePe Colour™.....	21
TePe Colour™ Compact.....	21
Brush Head Cap.....	21
TePe PlaqSearch™.....	22
Imprinting Service.....	23



**24–27**



**28–35**



**36–37**

## Special Products

---

TePe Universal Care™	25
TePe® Implant Orthodontic Brush	25
TePe® Bridge & Implant Floss	25
TePe Compact Tuft™	26
TePe Interspace™	26
TePe Special Care™	26
TePe Gentle Care™	26
TePe GOOD™ Tongue Cleaner	27
TePe® Denture Care	27
TePe Extra Grip	27

## Areas of Use

---

Healthy Gums	28
Implant Care	30
Orthodontic Care	32
Paedodontic Care	34

## Demo & Display

---

Mini Display	36
Interdental Brush Display	36
GOOD Table Display	36
Vanity Bag	36
Dental Model	36
Information Material	36
Medium Table Display	37
Table Display, Metal	37
Floor Stand, Metal	37
Product Presentation Case	37
Wall Display, Metal	37
Floor Stand, Cardboard	37


Product Guide, see pages 38–39.



# Interdental Cleaning

TePe has a long tradition of developing products for preventive dental care with special focus on interdental cleaning. Our range offers a variety of options to facilitate daily interdental self-care for all.





### TePe® Interdental Brush Original

The original range of TePe Interdental Brushes comprises nine colour coded sizes to fit narrow and wider interdental spaces. All brushes have plastic coated wire for safety and comfort. The user-friendly handle offers a stable grip which enables cleaning with controlled movements.

The smaller brush sizes (0–3) have a flexible neck for increased durability and access. The larger brush sizes (4–8) have a thicker bendable wire.



Interdental.




Missing tooth.



Furcation grade III.




Implants.




Orthodontic appliances.



Sizes 0–3: Curve the flexible neck for improved access between the posterior teeth.



Sizes 4–8: Curve the wire for improved access between the posterior teeth.





### TePe® Interdental Brush™ Extra Soft


This unique range of interdental brushes has extra soft filaments for gentle cleaning. It comprises six sizes, all with plastic coated wire for safe use. The colour coded sizes correspond to the original range, but in pastel colours.

Areas of use:

- After oral surgery
- Sore gingiva
- Dry mouth
- Sensitive teeth
- For everyone who prefers a softer brush



The plastic cap can be used for hygienic protection of the brush or to extend the handle.






## TePe Angle™

TePe Angle is designed for controlled and convenient interdental cleaning. It has an angled brush head for easier access, especially between the posterior teeth and from the palatal and lingual sides. The long and flat handle offers a stable and ergonomic grip. TePe Angle is a good choice for those who prefer a longer handle or when dexterity is reduced.


To fit both narrow and wider interdental spaces, TePe Angle is available in six colour coded sizes, corresponding to the original TePe Interdental Brush range. All sizes have plastic coated wire for gentle cleaning.



- Small, angled brush head for optimal reach
- Long, flat handle for a stable grip
- Plastic coated wire for safe use
- Six colour coded sizes




Posterior cleaning from the buccal side.




Posterior cleaning from the lingual side.



Anterior cleaning from the lingual side.



Cleaning of implants from the buccal side.




Cleaning of implants from the palatal side.

## Quality in every detail


Using an interdental brush is a highly effective way to remove dental plaque between the teeth. TePe's interdental brushes are developed in collaboration with dental experts to certify the highest efficiency, comfort and quality. They are continuously tested and evaluated to meet the demands of consumers and professionals worldwide.

The wide TePe Interdental Brush range offers an option for everybody – several sizes, two filament textures and two handle lengths with a straight or angled brush head. The sizes are colour coded for easy identification.

The brushes comply with the international standard for manual interdental brushes (ISO 16409) which stipulates nine sizes, 0-8. TePe offers brushes in all nine sizes to cover the whole span of interdental space varieties.



It is important to select the correct size. The brush must never be forced into a space.



All sizes and models have plastic coated wire for gentle cleaning.

Original									
Extra soft									
Angle									
Colour	Pink	Orange	Red	Blue	Yellow	Green	Purple	Grey	Black
ISO size (0–8)	0	1	2	3	4	5	6	7	8
Wire size (mm)	0.4	0.45	0.5	0.6	0.7	0.8	1.1	1.3	1.5
PHD (mm)	0.6	0.7	0.9	1.1	1.3	1.6	2.0	2.6	3.4

The ISO size is based on the measurement of the smallest hole which the brush can pass through, the passage hole diameter (PHD).

Interdental brushes are shown to be effective devices for interdental cleaning, more effective than dental floss or wooden sticks. Interdental brushes also have a positive effect on plaque scores, bleeding scores and probing pocket depth.

### Sälzer S, et al.

Efficacy of interdental mechanical plaque control in managing gingivitis – a meta-review.

J Clin Periodontol. 2015 Apr;42 Suppl 16:92-105.

# Accessories



## 1. Practice Box

Ideal for keeping TePe's interdental brushes within easy reach at the dental practice. Useful when trying out the correct size and instructing the patient. It can be filled with brushes from the multi-pack. A useful information chart with brush sizes is included.

The practice box for original and extra soft brushes has eight removable cases and the TePe Angle box has four cases. Each case holds 15 brushes.

## 2. Multi-pack

TePe's interdental brushes are available in single size packs of 25 brushes with individual plastic caps (grey and black without caps).

## 3. Travel Case

A practical travel case that fits six TePe Interdental Brushes original or extra soft. Airholes for good ventilation, allowing the filaments to dry.

## 4. Dental Model

The unique dental model has several educational features. It is useful for instructing cleaning techniques with TePe's interdental brushes and many other oral hygiene devices.

The model shows:

- Crowding
- Mesially inclined molar
- Missing tooth
- Furcation grade III
- Rotated, lingually inclined premolar
- Partially erupted third molar



**tepe interdental brushes**  
Developed in collaboration with dental experts

Original	Orange	Red	Yellow	Green	Purple	Black
1	2	3	4	5	6	7
0.5	0.6	0.7	0.8	0.9	1.0	1.1
1	2	3	4	5	6	7
0.5	0.6	0.7	0.8	0.9	1.0	1.1

1

2

## TePe® Gingival Gel with Chlorhexidine and Fluoride

This anti-bacterial gel combines the advantages of chlorhexidine (0.2%) and fluoride (0.32% NaF). Chlorhexidine has an antiplaque effect while fluoride prevents caries and protects sensitive root surfaces.

TePe Gingival Gel contains no abrasives, no alcohol and no colourants. The mint flavour and smooth formula enable pleasant and easy use with an interdental brush.


Areas of use:

- Antibacterial
- Antiplaque
- Caries protection

Thanks to its efficacy, chlorhexidine is considered the gold standard in the field of antiplaque agents.

**Varoni E, et al.**


Chlorhexidine (CHX) in dentistry: state of the art.  
Minerva Stomatol. 2012 Sep;61(9):399-419.



Drip gel along the interdental brush...



... or use the cap as a bowl to apply gel.



The gel stays between the teeth when the brush has been removed.




## TePe® Dental Tape

The dental tape has a flat, wide design for thorough plaque removal. It is coated with natural beeswax, which makes it glide easily between the teeth without shredding. Freshly mint flavoured with natural ingredients. The tape is made from strong and durable PTFE (PFAA free).

- Flat and wide
- Coated with natural beeswax
- Fresh mint flavour





## TePe GOOD Mini Flosser™

Pre-loaded dental floss holder for easy and convenient use. Strong, thin floss (PFAA free) that fits even in very narrow interdental spaces. The cleverly designed bite plane helps to gently insert the floss. The handle fits comfortably in the hand, and will from autumn 2019 be made from bio-based plastic.

- Bite plane for safe and easy use
- Small and slender handle
- Easy access to tight interdental spaces




Insert the floss by gently biting onto the bite plane.



Gently move the floss along each tooth surface.

TePe GOOD MiniFlosser™ will replace TePe Mini Flosser™ in autumn 2019.





- Silicone coating —
- Long working length —
- Gentle on the gums —
- Flexible and durable —
- Non-slip grip —


## TePe EasyPick™

TePe EasyPick™ encourages interdental cleaning on-the-go. It is durable, flexible and gentle on the gums. TePe EasyPick™ comes in two conical sizes with a steady, non-slip grip.

The entire cleaning part is covered with unique silicone lamellae. The wide lamellae in combination with the long working length offer efficient cleaning.




Easy to insert between teeth.



Flexible for easier access between the posterior teeth.

- Easy interdental cleaning
- Flexible and durable
- Use on-the-go



TePe EasyPick™ is available in two conical sizes, XS/S and M/L, to fit tight as well as wider interdental spaces.

**Ideal on-the-go**



A handy pocket case is included.


# Toothbrushing

TePe's user-friendly toothbrushes are developed for efficient and gentle plaque removal. A variety of models and sizes facilitate daily oral care for adults and children.



# TePe GOOD™ Toothbrush

A new Select toothbrush with the same TePe quality, design and feeling as always, but made from bio-based plastic. Thanks to the use of renewable raw materials, sugar cane and castor oil, we now manage to recirculate 95 % of the carbon dioxide emissions during the product's life cycle. It features a clever design with a tapered head and rounded filaments for great cleaning.



## TePe GOOD™ Regular Soft

This premium quality toothbrush with tapered brush head and ergonomic handle suits everyone.



## TePe GOOD™ Compact Soft

Recommended for children, but also for adults who prefer a smaller brush head. User-friendly handle which suits both children and adults.



## TePe GOOD™ Mini Extra Soft

TePe Mini is suitable from the first tooth until 3 years of age. The handle provides a stable grip for the adult brushing the child's teeth, but fits just as well in a small hand.



## TePe Select™

This premium quality toothbrush with tapered brush head and ergonomic handle suits everyone. Extra Soft, Soft or Medium.




## TePe Select™ Compact

Recommended for children, but also for adults who prefer a smaller brush head. User-friendly handle which suits both children and adults. Extra Soft, Soft or Medium.



## TePe Mini™

TePe Mini is suitable from the first tooth until 3 years of age. The handle provides a stable grip for the adult brushing the child's teeth, but fits just as well in a small hand. Extra Soft or Soft.



Filaments with rounded ends for gentle cleaning.



All toothbrushes can be angled after heating for better access.



TePe Select™



TePe Select™ Compact




TePe Mini™



## TePe Supreme™

Soft filaments in two levels for improved access. The longer filaments clean further in between the teeth and along the gum line, while the shorter filaments efficiently clean the rest of the tooth surfaces. Tapered brush head and a stable, non-slip handle. The neck can be angled without heating. Soft.



## TePe Supreme™ Compact

A toothbrush with a compact tapered brush head, suitable for everyone who prefers a smaller brush for easier access. Excellent for cleaning around orthodontic appliances. It has the same two-level filament design and non-slip handle as TePe Supreme. Soft.



## TePe Nova™

The easy-reach tip gives maximum access to difficult to reach areas, especially around the posterior teeth. Ergonomic, non-slip handle and tapered brush head. The neck can be angled without heating. Pink: Extra Soft. Yellow: Soft. Blue: Medium.





### TePe Graphic™

TePe Select with illustrations and coloured filaments. Soft.



### TePe Colour™


TePe Select with coloured filaments. Soft.



From 3 years

### TePe Kids™

Select Compact with illustrations, TePe Kids, is recommended from 3 years and up. Extra Soft or Soft.



### TePe Colour™ Compact

Select Compact is available with colourful filaments. Extra Soft.



0-3 years

### TePe Mini™

TePe Mini is available with charming prints created by Swedish Isa Form. Recommended from the first tooth. Extra Soft.




### Brush Head Cap

The practical cap suits most TePe toothbrushes and special brushes.



# Information online

In-depth product facts, user tips and pedagogic videos for professionals and patients are available on [www.tepe.com](http://www.tepe.com).



## TePe PlaqSearch™

PlaqSearch is a two-tone disclosing product to visualise plaque, colouring new plaque red and old plaque blue. It clearly shows the areas that need extra care. PlaqSearch is used at the dental practice for patient instruction and as an educational tool for teaching proper cleaning techniques. It is also convenient to use at home for self-assessed control and improved oral hygiene routines.

PlaqSearch is available as chewing tablets for use at the dental practice and at home.

- Educational and easy method to visualise plaque
- Differentiates new and old plaque
- Helps to prevent caries and gingivitis



Red shows new plaque and blue shows plaque older than 12 hours.

# Imprinting Service

Design your own toothbrush with TePe's imprinting service – an excellent marketing opportunity for dental practices.



Logotypes, images and messages can be printed on our Select and Select Compact toothbrushes. Choose between a variety of different handle colours and print colours for amazing designs.

# Special Products

TePe offers a wide range of cleaning devices developed for individuals with specific oral hygiene needs. Close collaboration with dental professionals ensures excellent function and a high level of quality. All products are clinically tested and evaluated.






## TePe Universal Care™

The unique, angled neck of Universal Care enables easy access to the lingual and palatal surfaces of implants. A small brush head with long filaments facilitates reach under the prosthetic construction. Universal Care is also suitable for retroclined teeth, lingual orthodontic appliances and other difficult-to-reach areas.



Lingual and palatal implant surfaces.



## TePe® Implant Orthodontic Brush


A long, thin neck and an extra slim brush head make this brush ideal for cleaning the buccal sides of implants and around orthodontic appliances.



Buccal implant surfaces.




Orthodontic appliances.




Extra slim brush head for improved access.

## TePe® Bridge & Implant Floss


A spongy filament floss specially designed for easy and efficient cleaning of implants, bridges and orthodontic appliances. The floss is equipped with stiff plastic ends for easy insertion.




Implant-supported bridge.



Single implant.



Bridge.



Orthodontic appliances.





## TePe Compact Tuft™


Small, angled brush head with a dense, rounded tuft. Ideal for precision cleaning of areas that need special attention.



Erupting teeth.



Along the gum line.



Orthodontic appliances.



Single implant.



Implant-supported bridge.




Attachments for overdentures.



## TePe Interspace™

An angled brush with a tapered brush tip for optimal access in difficult-to-reach areas.



Orthodontic appliances.




Furcations.




Distal surface of the last tooth.



Single implant.



Choose angle by attaching the brush head to either side of the neck.



Yellow: soft  
Blue: medium



## TePe Special Care™

Special Care with 12,000 ultrasoft filaments is recommended after oral surgery and for patients with extremely sore and delicate oral tissues. Two brush head sizes: regular (red), compact (blue).

Areas of use:

- After oral surgery
- Sore oral tissues
- After trauma
- Radiotherapy treatment



## TePe Gentle Care™

Gentle Care has super soft filaments, somewhat firmer than Special Care. Regular brush head.

Areas of use:

- After oral surgery
- Sensitive oral tissues
- After trauma
- Dry mouth
- Sensitive root surfaces

## TePe GOOD™ Tongue Cleaner

Using a tongue cleaner to remove the bacteria helps to prevent bad breath (halitosis) and to improve oral health. The carefully tested design of TePe GOOD™ Tongue Cleaner corresponds to the natural shape of the tongue and makes it easier to reach far back in the mouth. The slim shape reduces the gag reflex. The tongue cleaner has three slightly raised cleaning blades for a triple-effect, which means that one stroke usually is enough to clean the tongue efficiently. The tongue cleaner is now made from bio-based plastic.

- Removes bacteria with one stroke
- Prevents bad breath
- Dishwasher safe



## TePe® Denture Care

TePe Denture Care has extra long, strong filaments to ensure good accessibility and efficient cleaning. The brush is designed to resemble a conventional toothbrush.



Efficient cleaning of dentures.



Handle can be bent backwards for improved grip and access.



## TePe® Extra Grip

TePe Extra Grip assists persons with reduced manual dexterity, and fits most TePe toothbrushes and special brushes. It provides a comfortable, stable grip and weighs only 30 grams.




Insert/remove the brush by tapping the handle against for example, a table.

TePe Extra Grip is developed in collaboration with qualified occupational therapists at the Center for Technical Aid, Department of Rheumatology and Department of Hand Surgery at Skåne University Hospital, Sweden.




# Healthy Gums

Thorough oral hygiene is essential for maintaining healthy gums and preventing periodontal diseases. Motivating the patient and recommending the appropriate cleaning devices are key factors for efficient self-care.




## TePe Interdental Brushes

The wide range of sizes and designs enables individually adapted cleaning.




## TePe EasyPick™

An easy-to-use cleaning option in two conical sizes for tight to big interdental spaces.



## Dental Floss

Floss is suitable when the interdental space is narrow. Choose between easy-to-use TePe GOOD™ Mini Flosser or the flat, wide and waxed TePe Dental Tape.



## TePe® Gingival Gel

Antibacterial gel with chlorhexidine and fluoride. Recommended for a limited period of time as a complement to daily toothbrushing and interdental cleaning.



## TePe PlaqSearch™

A two-tone disclosing product to visualise plaque, colouring new plaque red and old plaque blue.

Effective plaque control provides the basis for maintaining gingival health. Toothbrushing and interdental cleaning form the cornerstones of mechanical plaque control.

### Chapple ILC, et al.

Primary prevention of periodontitis: managing gingivitis.  
J Clin Periodontol. 2015 Apr;42 Suppl 16:71-76.

### Kotsakis GA, et al.

A network meta-analysis of interproximal oral hygiene methods in the reduction of clinical indices of inflammation.  
J Periodontol. 2018 May;89(5):558-570.



# Implant Care

Proper care of implants is vital for long-term success. Professional guidance on good oral hygiene habits and suitable products is of great importance to prevent the development of peri-implant mucositis and peri-implantitis.



Peri-implant inflammation is common among implant patients. Good oral hygiene is of great importance, since poor oral hygiene increases the risk of development and progression of peri-implant disease.

**Derks J, Tomasi C**


Peri-implant health and disease. A Systematic review of current epidemiology.

J Clin Periodontol. 2015 Apr;42 Suppl 16:158-171.

**Heitz-Mayfield LJA, Salvi GE**

Peri-implant mucositis.

J Clin Periodontol. 2018 Jun;45 Suppl 20:237-245.



**TePe Universal Care™**

The angled neck facilitates cleaning from the lingual and palatal sides.




**TePe Implant Orthodontic Brush™**

Extra slim brush head, suitable for buccal cleaning of an implant-supported bridge.




**TePe Interdental Brushes**

Useful for inter-implant cleaning from the buccal and palatal sides, and under the connecting bar for overdentures.




**TePe EasyPick™**

Strong and flexible with silicone coating for efficient and gentle cleaning.



**TePe Compact Tuft™**

The small tuft facilitates cleaning of implants, and is also suitable around attachments for overdentures.



**TePe Bridge & Implant Floss™**

Spongy filament floss for cleaning around implants.



# Orthodontic Care

Fixed orthodontic appliances make cleaning more complicated. Plaque rapidly accumulates around brackets, bands and wires. To prevent gingivitis and caries, systematic cleaning is essential.



Fixed orthodontic appliances increase the number of sites where plaque can accumulate. Encouraging patients to maintain a high standard of plaque control is essential to reduce risk of caries and gingivitis.

**Sundararaj D, et al.**


Critical evaluation of incidence and prevalence of white spot lesions during fixed orthodontic appliance treatment: A meta-analysis.

J Int Soc Prev Community Dent. 2015 Nov-Dec;5(6):433-439.

**Srivastava K, et al.**


Risk factors and management of white spot lesions in orthodontics.

J Orthod Sci. 2013 Apr;2(2):43-49.



**TePe Supreme™**

The toothbrush with filaments in two levels facilitates access around orthodontic appliances.




**TePe® Implant Orthodontic Brush**

The extra slim brush head enables cleaning around brackets and wires.





**TePe Interspace™**

The tapered brush tip is ideal for cleaning around brackets and under wires.





**TePe Compact Tuft™**

The rounded tuft is suitable for cleaning around brackets and along the gum line.




**TePe Interdental Brushes**

Suitable for cleaning under the wire and for interdental cleaning under the retainer.



**TePe EasyPick™**

Easy interdental cleaning under the wire or the retainer.



**TePe® Bridge & Implant Floss™**

Efficient cleaning around brackets.



**TePe Universal Care™**

The angled neck facilitates cleaning of lingual orthodontic appliances.

# Paedodontic Care

Establishing good oral hygiene routines at a young age is the first step to a healthy mouth throughout life. Once the first tooth emerges, it is time to start using a toothbrush.





### TePe Mini™

A toothbrush with an extra small brush head, recommended from the very first tooth up to 3 years of age.

### TePe Select™ Compact

This toothbrush has a slightly larger brush head than TePe Mini. Suits children from 3 years and up.

### TePe Select™


The larger brush head makes Select suitable for older children.

### TePe Compact Tuft™


This single tufted brush is ideal for use on the critical occlusal surface on erupting molars.



Reaccumulation of plaque (coloured areas) 48 hours after professional toothcleaning.



Partially erupted molar.



Erupted molar with full chewing function.

Modified by TePe with permission from Prof. P. Axelsson, DDS, Odont. Dr.

### TePe PlaqSearch™

Two-tone disclosing tablet, indicating new and old plaque. Convenient tool for teaching proper brushing to children.

### TePe GOOD Mini Flosser™

Children with increased risk of caries and/or gingivitis need help to clean interdental. The Mini Flosser offers a stable grip for the adult helping the child.

### TePe Special Care™ Compact

A toothbrush with ultrasoft filaments and a small brush head. Recommended after trauma and oral surgery.

The quality of tooth brushing in children is dependent on parent support as well as toothbrushing technique.

**Collett B, et al.**

Observed child and parent tooth brushing behaviors and child oral health. *Int J Paediatr Dent.* 2016 May;26(3):184-192.

**Muller-Bolla M, Courson F**

Toothbrushing methods to use in children: a systematic review. *Oral Health Prev Dent.* 2013;11(4):341-347.

# Demo & Display

TePe offers a range of material for demonstration and information to aid chairside oral hygiene advice. Our product racks are perfect for display and storage.



## Mini Display

Designed for six toothbrushes and nine interdental brushes. Clear acrylic material.  
Item no 931 031.



## Interdental Brush Display

Suitable for 15 interdental brushes.  
Clear acrylic material.  
Item no 931 035.



## GOOD Table Display

Cardboard table display suitable for 2 SRP's (shelf ready packaging) with toothbrushes.  
Height: 29.5 cm Width: 18 cm. Depth: 16 cm.  
Item no 931086.



## Vanity Bag

Can be filled with products and information for your patient. Ideal as a tool for continued oral self-care after treatment.  
Height: 14 cm. Width: 25 cm. Item no 227.



## Dental Model

Educational dental model useful for teaching cleaning techniques with TePe's products.  
The picture illustrates an example of the dental model.



## Information Material

TePe provides a selection of instructive material for patients as well as dental professionals.



### Medium Table Display

Table display in cardboard with 8 hangers.  
Height: 63.5 cm. Width: 36.5 cm. Depth: 16 cm.  
Item no 932033.



### Table Display, Metal

Table display with rotating function, 8 hangers.  
Designed for all TePe products with hangers.  
Height: 53 cm. Width: 20 cm. Depth: 34 cm.  
Item no 931055.



### Floor Stand, Metal

Floor stand on wheels with 21 hangers, 1 shelf  
and 2 folder displays. Height: 164.5 cm.  
Width: 50 cm. Depth: 37.5 cm.  
Item no 931054



### Product Presentation Case

Includes the major part of the TePe range and is  
excellent for product demonstration.  
Dental Model is not included.  
Item no 952054.



### Wall Display, Metal

Wall display with 12 hangers and a shelf.  
Height: 80 cm. Width: 50 cm. Depth: 24.5 cm.  
Item no 931053.



### Floor Stand, Cardboard

Floor Stand made of cardboard with 16  
hangers. Height: 147 cm. Width: 50 cm.  
Depth: 33 cm. Item no: 931060.

# Product Guide



## TePe Interdental Brushes, pages 8–9

Available in different packages with a variety of languages.

Product	Multi-pack 25 pcs
0.4 mm pink	132 410
0.45 mm orange	132 420
0.5 mm red	132 430
0.6 mm blue	132 440
0.7 mm yellow	112 450
0.8 mm green	112 460
1.1 mm purple	112 470
1.3 mm grey	112 480
1.5 mm black	112 490
0.45 mm orange extra soft	122 425
0.5 mm red extra soft	122 435
0.6 mm blue extra soft	122 445
0.7 mm yellow extra soft	122 455
0.8 mm green extra soft	122 465
1.1 mm purple extra soft	122 475

## Interdental Brush Accessories, page 12

Product	Item No.	Item No.
Travel Case	922 028	922 428 with hanger
Dental Model	952050	
Practice Box Interdental Brushes	922 045 empty	922 045K filled
Practice Box Angle	922 041 empty	922 044 filled



## TePe Angle™, page 10

Product	Blister 6 pcs	Multi-pack 25 pcs
0.4 mm pink	154 610	154 410
0.45 mm orange	154 620	154 420
0.5 mm red	154 630	154 430
0.6 mm blue	154 640	154 440
0.7 mm yellow	154 650	154 450
0.8 mm green	154 660	154 460
Mixed pack Angle	154 699	




## Interdental Products, pages 14–17

Product	Item No.	Item No.
TePe EasyPick™, orange XS/S	242 310	
TePe EasyPick™, turquoise M/L	242 350	
TePe Gingival Gel™		982 373 blister
TePe Dental Tape™, 40 m	612 330	
TePe Mini Flosser™	621 100	
<b>NEW</b> TePe GOOD Mini Flosser™	622 153	



## Toothbrushes, pages 19–21

Product	Blister pack	Plastic bag
TePe Supreme™	352 688	352 100
TePe Supreme™ Compact	392 688	392 120
TePe Nova™, Extra Soft	312 687	312 110
TePe Nova™, Soft	312 688	312 120
TePe Nova™, Medium	312 689	312 130
<b>NEW</b> TePe GOOD™ Regular Soft	302 684	301 184
<b>NEW</b> TePe GOOD™ Compact Soft	302 685	301 185
<b>NEW</b> TePe GOOD™ Mini Extra Soft	302 686	301 186
TePe Select™, Extra Soft	322 687	322 110
TePe Select™, Soft	322 688	322 120
TePe Select™, Medium	322 689	322 130
TePe Colour™, Soft, pink bristles		322 183
TePe Colour™, Soft, green bristles		322 184
TePe Colour™, Soft, turquoise bristles		322 185
TePe Select™ Compact, Extra Soft	332 687	332 110
TePe Select™ Compact, Soft	332 688	332 120
TePe Select™ Compact, Medium	332 689	332 130
TePe Colour™ Compact, Extra Soft		332 182
TePe Mini™, Extra Soft	382 687	382 110



## Toothbrushes with illustrations, pages 20–21

Product	Blister pack	Plastic bag
TePe Graphic™, Soft	329 690	
TePe Kids™, Extra Soft	339 691	339 160
TePe Kids™, Soft	339 692	339 161
TePe Mini™, Extra Soft	382 687	
Brush Head Cap	922 024 <sup>box</sup>	

## Disclosing, page 22



Product	Tablets
TePe PlaqSearch™	992 676



## Special Products, pages 25–26

Product	Blister pack	Plastic bag
TePe Universal Care™	472 368	
TePe® Implant Orthodontic Brush	442 368	442 100
TePe® Bridge & Implant Floss	672 250 <sup>box</sup>	
TePe Compact Tuft™	432 368	432 100
TePe Interspace™, Soft	412 364	
TePe Interspace™, Medium	412 366	



## Special Products, pages 26–27

Product	Blister pack	Plastic bag
TePe Special Care™	452 357	452 150
TePe Special Care™ Compact		452 151
TePe Gentle Care™	462 368	462 150
<b>NEW</b> TePe GOOD™ Tongue Cleaner	942 385	
TePe® Denture Care	422 368	
TePe Extra Grip	912 210	



Made in  
Sweden

## We care for healthy smiles

TePe Munhygienprodukter AB is a Swedish company founded in 1965. The TePe range of interdental brushes, toothbrushes and dental sticks is developed in collaboration with dental experts with the aim of promoting good oral health. All design, development and production take place at the headquarters in Malmö. The products are marketed by subsidiaries in eight countries and distributors in over 60 countries, and available primarily through dental practices, pharmacies and retail. TePe has 340 employees and turnover of MSEK 750 (2018).

All trademarks herein are the property of TePe Munhygienprodukter AB.

### TePe Munhygienprodukter AB

Bronsåldersgatan 5  
213 76 Malmö, Sweden

Tel +46 (0)40 670 11 00

Fax +46 (0)40 670 11 11

[info@tepe.com](mailto:info@tepe.com)

[www.tepe.com](http://www.tepe.com)



Intertek

TePe is certified according to the international quality standard ISO 9001:2008 and the environmental standard ISO 14001:2004.

